

Product Code:

AL-E015

## CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

### UJT CHARACTERISTICS TRAINER

**OBJECTIVE:** To plot Static & Transfer Characteristics.

**FEATURES:**

Instrument comprises of Two DC IC Regulated Power Supplies 0-15VDC/150mA & 0- 30VDC/150mA.

Three MR65 Square/Digital meters for voltage & current measurement.

One UJT 2N2646 mounted behind the panel.

Instrument Fitted in Highly Shockproof Plastics Box.

Instrument Dimensions: 12"x8"

**ACCESSORIES:**

Sufficient no of 4mm Patch Cords for Interconnection

Instruction Manual for Student Workbook.

Power Cords 230v, 50Hz

**PHOTOGRAPH:**

