

Product Code:

AL-E004

CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

FET CHARACTERISTICS TRAINER

OBJECTIVE: To plot Static & Transfer Characteristics.

FEATURES:

Instrument comprises of Two DC IC Regulated Power Supplies 0-15VDC/150mA & 0- 3VDC/150mA.

Three MR65 Square/Digital meters for voltage & current measurement.

One FET BFW10 mounted behind the panel.

Instrument Fitted in Highly Shockproof Plastics Box.

Instrument Dimensions: 12"x8"

ACCESSORIES:

Sufficient no of 4mm Patch Cords for Interconnection

Instruction Manual for Student Workbook.

Power Cords 230v, 50Hz

PHOTOGRAPH:



AMBEYLABS