

Product Code:

AL-E013

CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

SCR CHARACTERISTICS TRAINER

OBJECTIVE: To plot Static & Transfer Characteristics.

FEATURES:

Instrument comprises of Two DC IC Regulated Power Supplies 0-5VDC/150mA & 0- 30VDC/150mA.

Three MR65 Square/Digital meters for voltage & current measurement.

One TYN604 mounted behind the panel.

Instrument Fitted in Highly Shockproof Plastics Box.

Instrument Dimensions: 12"x8"

ACCESSORIES:

Sufficient no of 4mm Patch Cords for Interconnection

Instruction Manual for Student Workbook.

Power Cords 230v, 50Hz

PHOTOGRAPH:

