

Product Code:

AL-E157

## CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

### THERMOCOUPLE TEMPERATURE TRANSDUCER TRAINER

#### SPECIFICATION :

Study of THERMOCOUPLE as temperature measuring transducer

THERMOCOUPLE probe with protection cover

3.5 digit DPM as temperature indicator

Bridge circuit for transducer

Differential amplifier with feedback

Lead compensation system

O/P provided on test points for monitoring & controlling.

Electrical Kettle & Glass thermometer for reference

Detailed instruction manual

#### PHOTOGRAPH:

