

Product Code:

AL-E154

CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

RTD TEMPERATURE TRANSDUCER TRAINER

SPECIFICATION :

Study of RTD as temperature measuring transducer
RTD (PT - 100) probe with protection cover
3.5 digit DPM as temperature indicator
Bridge circuit for transducer
Differential amplifier with feedback
Lead compensation system
O/P provided on test points for monitoring & controlling.
Electrical Kettle & Glass thermometer for reference
Detailed instruction manual

PHOTOGRAPH:

