

Product Code:

AL-E011

CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

DC REGULATED POWER SUPPLY 0-30V,1A/2A/3A

AMBEYLAB Digital Power Supplies has been designed and manufactured as per exact engineering standards using best quality components, offering overload/short circuit protection.

- Different controls for voltage and current controls
- Fine and coarse controls for variation
- Line Regulation is 0.01%. +10% Line variation
- & load regulation is 0.05% from no load to full load
- Our 3½ digit Digital Panel Meters are
- Compact panel meters for digital display of electrical quantities.
- The output is short circuit protected with crossbar electronics fold back current limiting circuit.

Applications:

These power supplies are very useful for R&D laboratories, control panel manufacturers, defense Establishment's school & college laboratories.

Technical specifications:

Output Voltage Output Current Meters
0-30V, 1A/2A/3A 2 DPM (3½ digit)

Output voltage:

0-30V $\pm 1\%$ continuously adjustable

Output current:

0-3A $\pm 1\%$ continuously adjustable

Operative voltage 230V, 50Hz $\pm 10\%$

Supplied with instruction manual

Picture:



Features

1. Trainer operated on 230V, 50 Hz \pm 10%
2. Supplied with instruction manual
3. Highly Short Circuit Proof.
4. With Digital Readout.