

Product Code:

AL-E002

## CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

### TRANSISTOR CHARACTERISTICS TRAINER

- Trainer used to study characteristics of Transistor of NPN/PNP in Common Emitter, Common Base,
- Kit comprises of 4 analog meters voltmeter 1V, 10V and ammeter of 250uA/50mA. 50mA of size 65mm with rectangular
- Built in regulated power supplies 0-1V & 10V required characteristics in different mode.
- Meter ranges are selectable with toggle switches
- 2 transistor NPN/PNP connected and their terminals are brought on front panel for easy connections via 4mm sockets.
- Front plate built with high class insulated Bakelite sheet with well printed circuits and symbols.
- Supplied with instruction manual.
- Operated on Mains power 230V, 50Hz.
- Supplied with sufficient nos. of patch cords & power inlet cable
- Supply provided through the 3 pin computer socket with appropriate fuse
- The trainer is housed in plastic cabinet to avoid shock
- Size of the trainer set 15"x8"

**Picture:**

