

Product Code:
AL-E311

CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

1. LOAD TEST ON 3 PHASE INDUCTION MOTOR:

Specifications:

- Electrical control panel with 1 Nos. 3 phase 3KVA transformer.(1KVA per phase)
- Panel size 24"X24".
- Electrical panel fitted with front panel of high class insulated Bakelite sheet 6mm of thickness.
- It is fitted with appropriate Digital voltmeter, ammeter to record the readings.
- Short circuit protection with the help of MCB.
- On board 3 Phase DOL Starter.
- Connections through 4mm highly shockproof safety sockets
- Supplied with instruction manual
- Digital Voltmeter Ammeter and Wattmeter are connected on the front panel to measurement the electrical voltage current system.

Product Photographs:

Experiments:

1. Load Test on 3 Phase Induction Motor
2. Speed Control on Three Phase Induction Motor.

Accessories:

- 1. 3 Phase Auto Transformer Close Type 5A per Phase**
- 2. 3 Phase Induction Motor 1Hp With Break Load Arrangement.**
- 3. Adequate No of 4mm Patch Cords for Interconnection.**
- 4. Instruction Manual**

