

Product Code:

AL-E419

## CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

### TWO PORT LADDER TRAINER

#### FEATURES:

Instrument comprises of DC IC Regulated Power Supplies 0-20V/150mA.  
On Board AC Supply 15V  
Two MR65 Square/Digital/Round meters for voltage & current measurement.  
T and PIE Network Provided on the board.  
Instrument Fitted in Highly Shockproof Plastics Box.  
Instrument Dimensions: 10"x8"

#### ACCESSORIES:

Sufficient no of 4mm Patch Cords for Interconnection  
Instruction Manual for Student Workbook.  
Power Cords 230v, 50Hz

#### PHOTOGRAPH:

