

**Product Code:**

AL-E156

## CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

### EXTRA HIGH TENSION (E.H.T.) POWER SUPPLY 10KV

**OUTPUT VOLTAGES:**

1. 0 to 5KV
  2. 0 to -5KV
  3. 0 to 10KV
  4. 6.3VAC/2A,
- High voltage proof up to 10KV via 4mm safety sockets

**LOAD CAPACITY:**

1. Max. 2mA (Short Circuit Current)
2. Max. 100uA ( Short Circuit Current)
3. Max. 200uA (Short Circuit Current)
4. 2A

**EXTERNAL VOLTAGES:**

1. 0-5VDC
2. 0-5Vp up to 1Hz.
3. Voltage Indicator: 2and Half Digit LED, 12.5mm.
4. Connection Voltage: 230V,50/60Hz.
5. Protection: T 0.5 for 230V; T 1.0 for 115V.
6. Power Consumption: 50VA.

**PHOTOGRAPH:**



AMBEYLABS