

Product Code:

AL-E309

## CATALOGUE

#118, Karanpuri, Behind B.D.Floor, Ambala Cantt, Haryana (India),

### POWER MEASUREMENT BY TWO WATTMETER METHOD

**SPECIFICATIONS:**

Panel size 42inches x24 inches.

Electrical panel consist of MS powder coated housing with back side removable.

Electrical panel fitted with front panel of high class insulated Bakelite sheet 6mm thick.

It is fitted with appropriate digital voltmeter, ammeter & wattmeter to record the readings.

Short circuit protection with the help of MCB.

Connections through 4mm safety terminals.

Supplied with instruction Manual

**PHOTOGRAPH:**

